

# No-Reference Quality Assessment Metric Protocol

Wenhui Zhu

**Date: 2022/10/11**


This is an evaluation criterion for unsupervised retinal image noise reduction tasks. Due to the lack of pairs of clean and noisy images, we formulate a metric for the evaluation of No-reference quality.

The main assessment indicators come from three points: (1) Lesion structure; (2) Background Color; (3) Generating extra Structure.


## (1) Lesion Structure

In the denoise task, it is easy for the generative model to ignore the lesion structure. The denoised images lesion and the original images do not match. **The main assessment method is the consistency of lesion structure and loss of lesion information.**

Here are some examples for reference



Lesion structure changed after denoised




Keep the lesion structure




Lesion structure loss

## (2) Background Color

We consider that many retinopathies are related to the background color, So the denoise images need to keep



Background color changed after denoise



Background color does not have big changes after denoise

the main color does not have a big change.

### (3) Generating extra Structures

For noise reduction tasks without paired data, most use adversarial generative models to denoise low-quality images by learning the distribution of high-quality images. However, a major problem with generative models is that generate some unexpected regions and structures. **For example, additional vessels, and new retinal regions.**

Generate new vessel structures

Generate new retinal regions



Keep structures consistent, no extra structures are added

